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[1] Lorber, M., A pharmacokinetic model for estimating exposure of Americans to dioxin-like compounds in the past, present, and future, *Sci. Tot. Environ.* 288 (2002), 81-95.

[2] United States Environmental Protection Agency (USEPA), Exposure and Human Health Reassessment of 2,3,7,8-Tetrachlorodibenzo-p-Dioxin and Related Compounds, Washington, D.C.:National Center for Environmental Assessment, U.S. Environmental Protection Agency (2000) EPA/600/P-00/001Be.

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[4] Council Directive 96/23/EC concerning the performance of analytical methods and the interpretation of results, *Off. J. Europ. Communities* L221/8, 17.8.2002.

[5] Council Directive 2002/69/EC, July 26, 2002 laying down the sampling methods and the methods of analysis for the official control of dioxins and the determination of dioxin-like PCBs in foodstuffs, *Off. J. Europ. Communities* L209/5, 6.8.2002.

[6] Council Directive 2002/70/EG, July 26, 2002 establishing requirements for the determination of levels of dioxins and dioxin-like PCBs in feedingstuffs, *Off. J. Europ. Communities* L209/15, 6.8.2002.

[7] Method 1613 Rev.B, Tetra- through Octa-Chlorinated Dioxins and Furans by Isotope Dilution HRGC/HRMS, U.S. Environmental Protection Agency Office of Water Engineering and Analysis Division, Washington, Oct 1994.

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

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[www.thermo.com](http://www.thermo.com)

Thermo Fisher Scientific, Madison, WI  
USA is ISO Certified.

©2007 Thermo Fisher Scientific Inc. All rights reserved. Golden Gate is a trademark of Specac. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

An80112\_ChS 01/07B